

# Perfect PAGE III

Written by **Lillian Zhu and Alan Tomusiak**

Checked by ?

For the Stanford Student Space Initiative Biology Subteam

October 29, 2017

## 1 Procedure Purpose

This experiment was run to determine whether DNA amount has an effect on the clarity of bands in a gel.

## 2 Overview

We ran a single gel with serial dilutions of oligos and ladders, looking for clarity and single base-pair resolution.

## 3 Safety Information

1. **SYBR Gold** is a mutagen and can penetrate laboratory gloves in a relatively short period of time, please change your gloves in the event of contamination.
2. Working in a communal lab space is dangerous. Do not assume your fellow workers cleaned up sufficiently

## 4 Materials

1. 1 15% PAGE Urea gel
2. Loading Buffer (dye)
3. TBE Buffer
4. 40 ng/uL 15 bp oligo
5. 40 ng/uL 16 bp oligo
6. Custom ladder mix
7. UltraPure water

## 5 Procedure

1. Rinse wells with 15  $\mu$ L TBE buffer. Pre-run for 20 minutes. Get the DNA samples out of the freezer.

### 5.1 Dilutions

- (a) 15 and 16 bp
  - For 1 in 5, add 10  $\mu$ L original sample to 40  $\mu$ L water.
  - For 1 in 10, add 10  $\mu$ L original sample to 90  $\mu$ L water.
  - For 1 in 100, add 10  $\mu$ L 10x dilution to 90  $\mu$ L water.

## (b) Ladder

- For 1 in 10, add 1  $\mu\text{L}$  original sample to 9  $\mu\text{L}$  water.
  - For 1 in 100, add 1  $\mu\text{L}$  original sample to 99  $\mu\text{L}$  water.
  - For 1 in 1000, add 1  $\mu\text{L}$  original sample to 999  $\mu\text{L}$  water.
  - For 1 in 10000, add 1  $\mu\text{L}$  10x dilution to 999  $\mu\text{L}$  water.
2. Obtain a sizable piece of parafilm. Pipette 5  $\mu\text{L}$  of Gel Loading Buffer in a row of 14 droplets.
  3. For the 10 bp ladder column, add 0.5  $\mu\text{L}$  ladder stock and 4.5  $\mu\text{L}$  TBE Buffer.
  4. For the remaining droplets, add 5  $\mu\text{L}$  of the appropriate sample.
  5. As you go, pipette up and down to mix thoroughly.
  6. Load the gels (with 10  $\mu\text{L}$  sample in each well) when they are finished pre-running. Ensure pipette tip is fully in the well, and depress slowly and carefully. Work quickly to minimize diffusion.

Well number	Sample
1	10 bp Ladder (control)
2	Undiluted 15 bp (200 ng)
3	1:5 15 bp (40 ng)
4	1:10 15 bp (20 ng)
5	1:100 15 bp (2 ng)
6	Undiluted 16 bp (200 ng)
7	1:5 16 bp (40 ng)
8	1:10 16 bp (20 ng)
9	1:100 16 bp (2 ng)
10	1:10 ladder mix
11	1:100 ladder mix
12	1:1000 ladder mix
13	1:10000 ladder mix
14	Blank

Figure 1: Wells

7. Run the gel at 150 V until dark blue dye is about 1 cm from the bottom. This should take about 1 hour.
8. Meanwhile, prepare 50 mL 1x SYBR Gold solution in an appropriate container.
9. After gel finishes, soak in SYBR for 30 minutes (try and be exact, but can probably vary within 5 min).
10. Proceed to imaging. Carefully position the gel, ensuring correct orientation, and optimize for intense bands.
11. Post resulting pictures and interpretations to Slack.