

The Long Weekend: Day 5 Review Analysis

Written by **Michael Uttmark**

Checked by —

For the Stanford Student Space Initiative Biology Sub-team

October 29, 2017

1 Procedure Purpose

To determine if reagents planned for use in a possible MVP (Minimally Viable Product) launch are viable after being left in solution for an extended period of time.

2 Overview

The reagents will be removed from various storage temperatures (room temp 23°C, -4°C and -20°C) and check for predicted florescence. They then will be combined with their “missing” components needed in order to preform a full Terminal Deoxynucleotidyl Transferase homopolymer extension. At this point, 6µL will be removed, heat denatured at 95°C, and stored at -20°C for future use. The samples will then be incubated and then heat denatured. The samples will then be examined via gel electrophoresis.

3 Safety Information

1. **SYBR Green I** is a mutagen and can penetrate laboratory gloves in a relatively short period of time, please change your gloves in the event of contamination. See <http://www.sigmaaldrich.com/MSDS/MSDS/DisplayMSDSPage.do?country=US&language=en&productNumber=S9430&brand=SIAL> for more information on the specifics of SYBR Green I.
2. **Ethidium Bromide** is a **serious mutagen** and is **significantly carcinogenic**. If working with considerable amounts, a **fume hood and respirator** are warranted. For more information see <https://www.sciencelab.com/msds.php?msdsId=9927667>
3. Working in a communal lab space is dangerous. Do not assume your fellow workers cleaned up sufficiently

4 Procedure

1. Remove the control samples from the room temperature, fridge and freezer storage.
2. Relabel samples if labels are no longer clear
3. Pipette samples into 96 well plate(s) in the order displayed below in Figure ??.
4. Run fluorescence tests using well plate viewer (Gain=50)
5. Pipette samples back into PCR tube

OPTIONAL STOP POINT

6. Pipette the following reagents into PCR tubes as described in Figure ??

	1	2	3	4	5	6	7	8	9	10	11	12
A	A1	H ₂ O	A2	H ₂ O	A3	H ₂ O
B	B1	H ₂ O	B2	H ₂ O	B3	H ₂ O
C	C1	H ₂ O	C2	H ₂ O	C3	H ₂ O
D	D1	H ₂ O	D2	H ₂ O	D3	H ₂ O
E	E1	H ₂ O	E2	H ₂ O	E3	H ₂ O
F	F1	H ₂ O	F2	H ₂ O	F3	H ₂ O
G	G1	J1	G2	J2	G3	J3
H	H1	I1	H2	I2	H3	I3

Figure 1: 96 Well Plate Sample Layout
(1: Room Temp, 2: -4°C, 3: 20°C)

	TdT	Sybr	Primer 1	Primer 2	5X TdT Buffer	dTTPs	Water	Ethidium Bromide
A	×	×	×	×	×	1μL	15.5μL	N/A
B	×	×	1.25μL *	1.25μL *	×	1μL	15.5μL	N/A
C	×	1μL *	1.25μL *	1.25μL *	×	1μL	15.5μL	N/A
D	×	1μL *	1.25μL *	1.25μL *	×	1μL	15.5μL	N/A
E	1μL	×	1.25μL *	1.25μL *	5μL	×	15.5μL	N/A
F	1μL	×	×	×	5μL	×	15.5μL	N/A
G	×	×	×	×	×	×	15.5μL	N/A
H	×	×	×	×	×	×	15.5μL	N/A
I	1μL	1μL *	1.25μL *	1.25μL *	5μL	1μL	×	N/A
J	1μL	N/A	1.25μL *	1.25μL *	5μL	×	15.5μL	×

Figure 2: Pipetting Instructions for PCR tubes

Note: * indicates the need for a dilution from stock, “x” indicates a reagent already present. See following instructions for dilutions.

Note: The dNTPs were originally supposed to be dTTPs, but due to human error, they were not. This is reflected in Figure ??.

Dilutions

4.1 SYBR Green I

- Vortex SYBR Green I stock
- Pipette 1μL SYBR Green I into a PCR Tube
- Pipette 9μL Water into solution
- Vortex until mixed

4.2 Primer 1

- Vortex Primer 1 stock
- Pipette 2.5μL Primer 1 into a PCR Tube
- Pipette 22.5μL Water into solution
- Vortex until mixed

4.3 Primer 2

- Vortex Primer 2 stock
- Pipette 2.5μL Primer 2 into a PCR Tube
- Pipette 22.5μL Water into solution
- Vortex until mixed

7. Incubate samples at 35°C for 60min
8. Stop the reaction by bringing the solution to 95°C for 10 minutes.
9. Then, bring to 2°C for 1 min for handling

OPTIONAL STOP POINT

10. Pipette samples A-I from all groups into 96 well plate(s)
11. Run florescence tests using well plate viewer

OPTIONAL STOP POINT

12. Prepare Agarose (Use Cobalt Chloride as a buffer for Gel)
 - (a) To volume 80mL Buffer, add 2.2g Agarose (2.75% Gel).
 - (b) Add 8µL of SYBR Green to Agarose once lukewarm
 - (c) Refill Uytengsu Buffer stock
13. Prepare Primer Mix
 - (a) Pipette 1µL Primer 1 to PCR tube
 - (b) Pipette 1µL Primer 2 to solution
 - (c) Pipette 78µL Water to solution
 - (d) Vortex to Mix
14. Plan Sample location in gels (see tables below)
15. Add Ladder to Gel (1µL of loaded solution)
16. Add samples and Primer Mix to gel (5µL sample + 1µL Loading dye \Rightarrow 5µL of loaded solution)
17. Run Gel at 100V for 20 minutes
18. Analyze Gel on Gel Viewer
19. Run Gel at 100V for 60 more minutes for a total of 80 minutes
20. Analyze Gel on Gel Viewer
21. Post pictures to SSI Slack

Stop Procedure

1. Pipette samples into PCR tubes if not already contained in an appropriate manner
2. Label containers
3. Freeze samples at -20°C

Sample location, top to bottom (wells on the right)	Gel 1 (Room Temp)	Gel 2 (-4°C)	Gel 3 (-20°C)	Gel 4
1	G2	H1	H3	J3
2	G1	G3	H2	J2
3	B3	D3	F3	J1*
4	B2	D2	F2	I3
5	B1	D1	F1	I1
6	A3	C3	E3	I2
7	A2	C2	E2	
8	A1	C1	E1	
9	Primer	Primer	Primer	Primer
10	Ladder	Ladder	Ladder	Ladder

For a detailed description of the gels, see section ??: Gels & Accompanying Analysis.

Figure 3: Gel Table

5 Overall Analysis

These data gathered by this experiment support the conclusion that the conditions present in samples A-C may result in the degradation of DNA over relatively short periods of time (13 days) at room temperature, -4°C and -20°C. However, the conditions in samples A-C still seem to be conducive to Terminal Deoxynucleotidyl Transferase and its continued activity after short periods of time at room temperature, -4°C and -20°C. Further, there seems to be evidence that, while the results from samples D, F and G are faint (in all strata), the conditions present in samples D, F and G seem to be conducive to Terminal Deoxynucleotidyl Transferase and its continued activity after short periods of time at room temperature, -4°C and -20°C. Nothing can be said about samples E and H from any of the gels. Further, these data support the hypothesis that Terminal Deoxynucleotidyl Transferase works better when stored at -20°C rather than -4°C and better at -4°C than at room temperature (this relation is transitive). However, the quality control samples (the Is) indicate that there was a measurable amount of contamination from some source, probably human error. As such, while the results from A-C still seem to be indicative of truth (due to their strong signal and consistent trends), results derived from samples D, F and G are reasonably considered to be questionable and will need further verification in the future if their implications become paramount.

6 Florescence Analysis

	Sample	Room Temp.	-4°C	-20°C	Deviation	$[10^{-2}]$	Room Temp.	-4°C	-20°C
Green	A	5.07	-1.41	-2.77	4.19		0	0	0
	F	1657.87	1657.87	28	N/A		16	16	0
	G	1657.87	1657.87	529.48	N/A		16	16	5
	H	1657.87	1657.87	1657.87	N/A		16	16	17
Yellow	J	1.98	37.22	-0.72	21.17		0	0	0
Red	B	2.81	-2.8	-3.7	3.52		0	0	0
	C	83.16	39.46	9.05	37.25		0	0	0
	D	145.85	18.46	2.83	78.45		1	0	0
	E	771.85	340.93	18.18	378.13		7	3	0
	I	2.68	2.39	1.47	0.63		0	0	0

Figure 4: Aggregated data pertaining to the florescence of samples (see Section?? for more information). Overflows replaced with 1657.87, the highest recorded value.

6.1 SYBR Green I Related Sample Analysis

Samples A, F, G and H were all expected to fluoresce when excited with 485 nm light (close to 497 nm, which is SYBR Green I's peak absorbance frequency for fluorescence [?]). These samples were expected to fluoresce because they contained SYBR Green I and at least one Primer. Some samples (B and E) contained SYBR Green I or a Primer, but not both, and were therefore not expected to fluoresce. Accordingly, samples F, G and H fluoresced significantly at Room Temperature, with their fluorescence tapering off from Room Temperature to -4°C and from -4°C to -20°C. This indicates that the change in temperature is positively correlated with the fluorescence of SYBR Green I in each

samples respective storage conditions. However, it is important to note that F1, F2, G1, G2, H1 and H2 all overflowed. This missing data might contradict this observed trend. More data is needed for strong support to be established. As this data stands, little can be said for sure.

Further, sample A did not fluoresce at all, and had a fluorescence pattern statistically identical to that of water. Either the sample was prepared incorrectly or this sample points to a severe instability of SYBR Green I in sample A's conditions. More data is needed for strong indication to be established and as this data stands, little can be said with any conviction. Still, this indicates that another experiment pertaining to A's conditions may be meritorious to complete. You can see this trend in Figure ??.




Figure 5: Storage Temperature vs. Normalized Arbitrary Florescence Units

Additionally, sample E fluoresced significantly even though it was supposed to contain DNA. This indicates one of two things: (1) there was an issue with contamination or (2) the instability of SYBR Green I caused this peculiar behavior. With the data at hand, it seems much more likely that there was an issue with contamination instead of an odd facet of SYBR Green's instability marvelously producing this behavior. With only three data points, little can be said for sure. D's measured florescence was insignificant at 1%.

6.2 Ethidium Bromide Related Sample Analysis

The data in Figures ?? and ?? indicates that only J2 had active Ethidium Bromide. Moreover, only J2 fluoresced significantly when exposed to 485 nm light during the SYBR Green I assays. There are multiple things that could explain this, but most fall flat:

1. Human Error: We did not add Ethidium Bromide to J1 or J3. This seems unlikely, as Ethidium Bromide has a very saturated color in solution, making it hard to pass over.
2. Ethidium Bromide degraded over time due to the temperature. This is unlikely, as J1 and J3 are the samples that didn't fluoresce, which implies that temperature isn't a particularly good predictor of fluoresce (however, one must always remember the imperative difference between correlation and causation).
3. The conditions of Ethidium Bromide caused degradation over time. This is also unlikely, as all three contained the same reagents.
4. Due to a human error, some reagent that was added to J2 was not added to J1 and/or J3 resulting in J2's anomalous behavior. This conclusion seems most likely, but more work should be done to support this hypothesis beyond our current data set with a cardinality of three.

6.2.1 Ethidium Bromide Sample Assay

Regrettably, after the procedure was completed and samples disposed of, it was discovered that the Ethidium Bromide samples were examined via absorbance, not florescence. Moreover, the samples were examined at 258 nm rather than 360 nm which is the peak absorbency frequency for Ethidium Bromide [?]. However, as we can see from Figure ??, there is still a significant absorbance of light at 258 nm (roughly 40%). As such, we can still determine that the Ethidium Bromide should still be detectable via absorption data. Moreover, it's important to note that the solvent (water) does not have a significant absorbance at 258 nm [?] and therefore should not impact our results significantly (See Figure ??).




Figure 6: Measured ultraviolet absorption spectra of liquid water[?]




Figure 7: Absorption spectrum of Ethidium Bromide in water at room temperature

Image by Mark Somoza (Own work) [GFDL, CC-BY-SA-3.0 or CC BY 2.5], via Wikimedia Commons

7 Gels & Accompanying Analysis

Continued on next page.




Figure 8: Gel 1

7.1 Gel 1 Analysis

As one can clearly see above (or rather not see) samples A and B are not visible. Further, there is no residue in their wells (normal activity for most DNA samples). These two facts indicate that these samples either were not added to the gel or did not contain DNA. (It is extremely unlikely that SYBR I would have stained the other samples and excluded A and B). We see extension in sample G2, but not in G1. This indicates that G1 did not see extension at room temperature, even though it was had all the needed reagents except for water. Further, this provides evidence that storing a "full pot" of reagents at room temperature leads to inactive reagents after 5 days, which is a relatively short period of time. However, this also provides support for the conclusion that storage at -4 degrees Celsius rather than at room temperature can in fact make a difference in reagent performance. (More images of this gel below.)

	1	2	3	4	5	6	7	8	9	10	11	12
A	585.5	554	459.5	484	433	440.5	417.5	447.5	430	449	454	448.5
B	541.5	545	432.5	449.5	415	434.5	428.5	411	438.5	433	442.5	459.5
C	2104.5	460.5	1254.5	445	663	429.5	434.5	409.5	435	421.5	454.5	450
D	3324	457	846	468	542	439.5	412	438.5	420.5	440	443	452
E	15500.5	461	7118.5	439	840.5	427	421.5	461.5	419.5	435	465.5	464.5
F	OVRFLW	509	OVRFLW	496.5	1031.5	451	435.5	441	425.5	435	477	451.5
G	OVRFLW	1.693	OVRFLW	3.724	10786	1.679	454	454	437	467.5	477	476
H	OVRFLW	539	OVRFLW	533.5	32735	515.5	462.5	458	467	468	493	464.5

Figure 9: Averages of the values from Figure ?? and Figure ?? with Ethidium Bromide sample related values (G2, G4 and G6) replaced with values from Figure ?. Please note that the Ethidium Bromide related values are in different units, see Section ?? for more information.

	1	2	3	4	5	6	7	8	9	10	11	12
A	589	559	451	480	450	437	401	449	423	444	457	443
B	555	538	435	429	412	452	422	410	445	456	441	457
C	2123	480	1297	450	663	436	449	418	445	417	468	465
D	4665	471	799	469	548	427	417	417	421	430	450	450
E	15583	467	13244	429	1013	433	414	472	414	450	466	463
F	OVRFLW	499	OVRFLW	470	1053	443	420	444	416	420	472	437
G	OVRFLW	504	OVRFLW	1234	10934	490	453	449	439	464	478	463
H	OVRFLW	548	OVRFLW	536	34587	541	463	460	469	461	497	444

Figure 10: Read One: Excitation: 485 nm, Emission: 528 nm (SYBR Green I Fluorescence Assay)

	1	2	3	4	5	6	7	8	9	10	11	12
A	582	549	468	488	416	444	434	446	437	454	451	454
B	528	552	430	470	418	417	435	412	432	410	444	462
C	2086	441	1212	440	663	423	420	401	425	426	441	435
D	1983	443	893	467	536	452	407	460	420	450	436	454
E	15418	455	993	449	668	421	429	451	425	420	465	466
F	OVRFLW	519	OVRFLW	523	1010	459	451	438	435	450	482	466
G	OVRFLW	547	OVRFLW	1188	10638	456	455	459	435	471	476	489
H	OVRFLW	530	OVRFLW	531	30883	490	462	456	465	475	489	485

Figure 11: Read Two: Excitation: 485 nm, Emission: 528 nm (SYBR Green I Fluorescence Assay)

	1	2	3	4	5	6	7	8	9	10	11	12
A	1.72	1.706	1.661	1.698	1.666	1.661	1.693	1.697	1.706	1.75	1.763	1.773
B	1.681	1.653	1.661	1.661	1.625	1.647	1.632	1.642	1.663	1.696	1.723	1.753
C	1.754	1.682	1.686	1.662	1.621	1.644	1.655	1.654	1.729	1.709	1.716	1.744
D	1.721	1.663	1.641	1.635	1.623	1.618	1.633	1.626	1.646	1.678	1.714	1.746
E	2.015	1.67	1.667	1.636	1.617	1.61	1.625	1.649	1.68	1.689	1.761	1.78
F	1.997	1.757	1.953	1.721	1.68	1.685	1.68	1.705	1.736	1.749	1.799	1.832
G	2.702	1.693	2.218	3.709	1.669	1.678	1.691	1.696	1.724	1.739	1.771	1.829
H	1.895	1.809	2.726	1.839	1.801	1.722	1.719	1.762	1.782	1.767	1.858	1.881

Figure 12: Read One: Absorbance at 285 nm (Ethidium Bromide Assay, see ?? for more information)

	1	2	3	4	5	6	7	8	9	10	11	12
A	1.719	1.706	1.662	1.697	1.664	1.659	1.692	1.697	1.705	1.751	1.764	1.773
B	1.68	1.654	1.661	1.66	1.623	1.646	1.631	1.642	1.663	1.695	1.723	1.752
C	1.757	1.683	1.685	1.661	1.621	1.643	1.655	1.654	1.729	1.708	1.715	1.743
D	1.719	1.663	1.639	1.634	1.622	1.618	1.633	1.628	1.649	1.679	1.714	1.745
E	2.016	1.669	1.667	1.636	1.616	1.609	1.626	1.649	1.68	1.688	1.76	1.78
F	1.987	1.757	1.939	1.721	1.68	1.686	1.681	1.707	1.734	1.749	1.801	1.833
G	2.667	1.693	2.129	3.724	1.668	1.679	1.691	1.697	1.723	1.738	1.771	1.831
H	1.879	1.808	2.613	1.841	1.793	1.722	1.72	1.763	1.781	1.768	1.858	1.882

Figure 13: Read Two: Absorbance at 285 nm (Ethidium Bromide Assay, see ?? for more information)